




SURGICAL
CATALOGUE
EDITION I

SCHAFTARTEN TYPES OF SHANK

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1456 Trephine Basic Kit 6

1457 Trephine Kit Long 7

1458 Biopsy Punch Set for Extraction of Tissue Grafts 8

Periimplantitis 8 - 10

- EN Periimplantitis
- FR Périimplantite
- ES Periimplantitis

1462 Carbide Finishers for Periimplantitis Treatment 8

1912 Periimplantitis/Paradontitis Kit 9

1454 Degranulation Kit 10

Alveoplastik 11

- EN Alveoplasty
- FR Alvéoplastie
- ES Alveoplastia

1471 Alveoplasty Surgical Kit 11

KITS AUSWAHL
KITS ASSORTMENT

1459

P. 4

1461

P. 4

1470

P. 5

1463

P. 5

1456

P. 6

1457

P. 7

1458

P. 8

1462

P. 8

1912

P. 9

1454

P. 10

1471

P. 11

SINUSLIFTING

SINUS LIFT

Kit 1459

LATERAL APPROACH KIT FOR EXTERNAL SINUS ELEVATION HP

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The Surgical kit 1459 contains instruments with handpiece shank for the lateral approach in an external sinus lift. Carbide cutters in three diameters and diamond instruments in three grains allow for the gradual thinning of the facial bone wall of the maxillary sinus down to the Schneiderian membrane. The large diameter of the instruments protects the membrane from unintentional perforation.

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The Surgical kit 1461 contains instruments with right-angled shank for the lateral approach in an external sinus lift. Carbide cutters in three diameters and diamond instruments in two grains allow for the gradual thinning of the facial bone wall of the maxillary sinus down to the Schneiderian membrane. The large diameter of the instruments protects the membrane from unintentional perforation.

SHAPE	FIGURE	SHANK	REF	SIZE	PIECES
	C141F	HP	500 104 001 251	018	1
	C141F	HP	500 104 001 251	031	1
	C141F	HP	500 104 001 251	050	1
	CX161R	HP	500 104 408 298	018	1
	801	HP	806 104 001 524	050	1
	801G	HP	806 104 001 534	050	1
	801SG	HP	806 104 001 544	050	1

SHAPE	FIGURE	SHANK	REF	SIZE	PIECES
	C141F	RAL	500 205 001 251	018	1
	C141F	RAL	500 205 001 251	031	1
	C141F	RAL	500 205 001 251	050	1
	CX161R	RAL	500 205 408 298	018	1
	801	RA	806 204 001 524	050	1
	801SG	RA	806 204 001 544	050	1

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The Surgical Kit 1470 contains a selection of conical Tungsten Carbide Cutters with handpiece shank for a broad field of application in oral surgery. The cutters have a particularly easy cutting and efficient saw toothing and a specifically designed twist angle. Further, the narrow working parts allow for an optimal view on the working area.

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The Bone Saw Kit provides a selection of diamond saws in the three most important diameters. They are used for the application in bone spreading or splitting operations. Due to the various sizes, the set offers the clinician the possibility of having a proper selection of these saws during surgery.

SHAPE	FIGURE	SHANK	REF	SIZE	PIECES
	CX162A	HP	500 104 409 019	016	1
	C163	HP	500 104 406 297	014	1
	C166	HP	500 104 409 297	021	1
	C254	HP	500 104 415 296	012	1
	C408M	HP	500 104 409 338	016	1
	C165	HP	500 104 408 297	023	1

SHAPE	FIGURE	SHANK	REF	SIZE	PIECES
	231DEF	RA	806 204 064 504	070	1
	231DEF	RA	806 204 064 504	100	2
	231DEF	RA	806 204 064 504	130	1

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The Trephine Basic Kit offers a selected range of high quality trephines in four diameters. These tools are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing implants. The trephine drills have a depth of 10 mm and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored. Furthermore, the pre-drilling trephines in the same diameter included in the set enable a safe application of the trephines.

SHAPE	FIGURE	SHANK	REF	SIZE	LENGTH MM			PIECES
	DV229	RA	330 205 557 001	040	2.0	4.0	4.7	1
	DV229	RA	330 205 557 001	050	2.0	5.0	5.7	1
	DV229	RA	330 205 557 001	060	2.0	6.0	6.7	1
	DV229	RA	330 205 557 001	070	2.0	7.0	7.7	1
	229RF	RAL	330 205 486 001	040	10.0	4.0	5.0	1
	229RF	RAL	330 205 486 001	050	10.0	5.0	6.0	1
	229RF	RAL	330 205 486 001	060	10.0	6.0	7.0	1
	229RF	RAL	330 205 486 001	070	10.0	7.0	8.0	1

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The Trephine Kit Long offers a selected range of high quality trephines in four diameters. These tools are designed specifically for creating cylindrical bone grafts quickly and accurately and for removing implants. The trephine drills have a depth of 14 mm* and have been developed specifically with mandibular surgery and implantology in mind. Thanks to the well-defined laser depth marking in 2 mm intervals, the drilling depth can accurately be monitored. Furthermore, the pre-drilling trephines in the same diameter included in the set enable a safe application of the trephines.

SHAPE	FIGURE	SHANK	REF	SIZE	LENGTH MM			PIECES
	DV229	RA	330 205 557 001	040	2.0	4.0	4.7	1
	DV229	RA	330 205 557 001	050	2.0	5.0	5.7	1
	DV229	RA	330 205 557 001	060	2.0	6.0	6.7	1
	DV229	RA	330 205 557 001	070	2.0	7.0	7.7	1
	229XLRF	RAL	330 205 556 001	040	18.0	4.0	5.0	1
	229XLRF	RAL	330 205 556 001	050	18.0	5.0	6.0	1
	229XLRF	RAL	330 205 556 001	060	18.0	6.0	7.0	1
	229XLRF	RAL	330 205 556 001	070	18.0	7.0	8.0	1

PERIIMPLANTITIS PERIIMPLANTITIS

Kit 1458

BIOPSY PUNCH SET FOR EXTRACTION OF TISSUE GRAFTS

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The Punch Basic Kit contains Biopsy Punches in five different diameters for the performance of tissue removal. They serve for the application in immediate implant placement, for taking mucosa grafts in the scope of soft tissue management or biopsy as well as for uncovering the gingiva after implantation. The Biopsy Punches create clean cut surfaces and care for minimal tissue traumatization.

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The Periimplantitis Kit contains perfectly matched tungsten carbide finishing burs in egg and flame shapes for the idea peri-implantitis treatment. All shapes come in standard and ultra-fine toothings. The shapes, sizes and toothings of the working parts are perfectly suited for the intraoral treatment of titanium and for different shapes of implant neck and shoulder. All finishers have an extra-long FG shank (total length 32 mm) which allows particularly the treatment of difficult-to-access, deep areas. For the peri-implantitis treatment, the finishers are used clockwise and are guided around the exposed implant counter-clockwise so the implant surface is smoothed.

SHAPE	FIGURE	SHANK	REF	SIZE	LENGTH MM	□	□	PIE- CES
	225RF	RA	330 204 485 373	030	6.0	3.0	3.7	1
	225RF	RA	330 204 485 373	050	6.0	5.0	5.7	1
	225RF	RA	330 204 485 373	060	6.0	6.0	6.7	1
	225RF	RA	330 204 485 373	070	6.0	7.0	7.7	1


SHAPE	FIGURE	SHANK	REF	SIZE	PIECES
	C379	FGXXL	500 317 277 072	014	1
	C379	FGXXL	500 317 277 072	023	1
	C48L	FGXXL	500 317 249 072	014	1
	C48L	FGXXL	500 317 249 072	023	1
	C379U	FGXXL	500 317 277 032	014	1
	C379U	FGXXL	500 317 277 032	023	1
	C48LU	FGXXL	500 317 249 032	014	1
	C48LU	FGXXL	500 317 249 032	023	1

Kit 1912

PERIIMPLANTITIS/PARADONTITIS

DEVELOPED WITH PD DR. SAHRMANN (BASEL, SWITZERLAND)

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
A professionally compiled assortment to ensure sustainable success in implantoplasty. The instruments, defined in collaboration with PD Dr. Sahrman, support controlled osteoplasty and smoothing and polishing of exposed titanium implant structures. With this set, exceptionally smooth surfaces are achieved, risks for progressive peri-implantitis are significantly reduced, and the preservation of implants is promoted.

SHAPE	FIGURE	SHANK	REF	SIZE	PIECES
	638.ARK	RA	635 204 110 505	025	1
	652.GRN	RA	655 204 173 523	035	1
	9501M	RA	653 204 243 513	030	3
	9831	RA	803 204 243 502	030	1


SHAPE	FIGURE	SHANK	REF	SIZE	PIECES
	811LG	FG	806 314 039 534	037	1
	C1S	RAL	500 205 001 003	016	1
	C1S	RAL	500 205 001 003	023	1

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
Osteoplastik mittels Rosenbohrer **C1S**
Osteoplasty with the round **C1S** bur

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
Degranulation kit is excellently suitable for further bone regeneration procedures. Degranulation burs can be used in implantology, periodontal surgery, endodontic surgery, and root removal. Very rough diamond edges allow fast removal of soft granulation tissue inaccessible to manual treatment. Diamond instruments do not damage bone structure compared to carbide burs; therefore safe when you work near the anatomical structures. The degranulation kit 1454 consists of four degranulation drills with diameter 1.0 mm, 2.5 mm, 3.0 mm, 3.5 mm.


SHAPE	FIGURE	SHANK	REF	SIZE	PIECES
	801LD	RAXL	806 206 697 534	010	1
	801LD	RAXL	806 206 697 534	025	1
	801LD	RAXL	806 206 697 534	030	1
	801LD	RAXL	806 206 697 534	035	1

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The Alveoplasty Surgical Kit contains the most important instruments for surgical contouring and modification of the alveolar ridge when unwanted unevenness of the jawbone occurs after tooth extraction or bone degeneration. In this case, the kit offers the best way to smooth the jawbone. The Alveoplasty Kit contains initial bur and two Lindemann burs, which allow for the performance of precise, initial osteotomies, even in hard cortical bone. To perform the alveoplasty itself, carbide cutters and specially designed diamond and ridge contourers are included. The rounded working parts optimally protect the soft tissue against injuries. Thus, the alveolar ridge is prepared carefully and precisely for implant insertion or denture restoration.

SHAPE	FIGURE	SHANK	REF	SIZE	PIECES
	186RF	RA	330 204 684 377	018	1
	C166	RAL	500 205 409 297	021	1
	165RF	RAL	330 205 408 297	023	1
	C71	HP	500 104 001 175	040	1
	CX77	HP	500 104 237 190	060	1

SHAPE	FIGURE	SHANK	REF	SIZE	PIECES
	801	RA	806 204 001 524	033	1
	D8307XG	HPS	806 103 307 554	060	1
	D8307XG	HP	806 104 307 554	090	1
	RCB00	RAL	-	075	1

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



2 2

APPLICATION

CHIRURGISCHE INSTRUMENTE AUS STAHL
SURGICAL INSTRUMENTS MADE OF STEEL

1RF

Kugelförmig (rund), Stahl rostfrei
spherical (round), stainless steel

Handstück | Straight Handpiece


FIGURE	SHANK	REF	Ø						
1RF	HP	330 104 001 001	007	009	010	012	014	018	023
			5	5	5	5	5	5	5
APPLICATION & HYGIENE									

141RF

6 Schneiden, Stahl rostfrei
6 blades, stainless steel

Handstück | Straight Handpiece



FIGURE	SHANK	REF	Ø				
141RF	HP	330 104 001 291	023	025	027	031	035
Winkelstück Right Angle							
141RF	RAXL	330 206 001 291	023				
			2	2	2	2	2
APPLICATION & HYGIENE							

162

Knochenfräser
bone cutter

Turbine | Friction Grip

FIGURE	SHANK	REF	Ø
162	FGXL	310 316 408 297	016
Winkelstück Right Angle			
162	RAL	310 205 408 297	016
		L mm	9,0
			2
APPLICATION & HYGIENE			

162RF

Knochenfräser, Stahl rostfrei
bone cutter, stainless steel

Handstück | Straight Handpiece

FIGURE	SHANK	REF	Ø
162RF	HP	330 104 408 297	016
Winkelstück Right Angle			
162RF	RAL	330 205 408 297	016
Turbine Friction Grip			
162RF	FGXL	330 316 408 297	016
		L mm	9,0
			2
APPLICATION & HYGIENE			

163RF

Knochenfräser, Stahl rostfrei
bone cutter, stainless steel


Handstück | Straight Handpiece

FIGURE	SHANK	REF	Ø
163RF	HP	330 104 406 297	014
		L mm	5,0
			2
APPLICATION & HYGIENE			

164RF

Knochenfräser, Stahl rostfrei
bone cutter, stainless steel


Handstück | Straight Handpiece

FIGURE	SHANK	REF	Ø
164RF	HP	330 104 407 297	018
		L mm	6,0
			2
APPLICATION & HYGIENE			

165RF

Knochenfräser, Stahl rostfrei
bone cutter, stainless steel

Handstück | Straight Handpiece

FIGURE	SHANK	REF	Ø
165RF	HP	330 104 408 297	023
Winkelstück Right Angle			
165RF	RAL	330 205 408 297	023
		L mm	7,0
			2
APPLICATION & HYGIENE			

166

Knochenfräser, konisch spitz
bone cutter, conical pointed

Handstück | Straight Handpiece

FIGURE	SHANK	REF	Ø
166	HP	310 104 409 297	021
Winkelstück Right Angle			
166	RAL	310 205 409 297	021
		L mm	10,0
			2
APPLICATION & HYGIENE			

166RF

Knochenfräser, konisch spitz, Stahl rostfrei
bone cutter, conical pointed, stainless steel

Handstück | Straight Handpiece

FIGURE	SHANK	REF	Ø
166RF	HP	330 104 409 297	021
		L mm	10,0
			2
APPLICATION & HYGIENE			

167RF

Knochenfräser, Stahl rostfrei
bone cutter, stainless steel

Handstück | Straight Handpiece

FIGURE	SHANK	REF	Ø
167RF	HP	330 104 410 297	023
		L mm	10,0
			2
APPLICATION & HYGIENE			

168RF

Knochenfräser, Stahl rostfrei
bone cutter, stainless steel

Handstück | Straight Handpiece

FIGURE	SHANK	REF	Ø
168RF	HPL	330 105 411 297	023
		L mm	22,0
			2
APPLICATION & HYGIENE			

168

Knochenfräser
bone cutter

Handstück | Straight Handpiece

FIGURE	SHANK	REF	Ø
168	HPL	310 105 411 297	023
		L mm	22,0
			2
APPLICATION & HYGIENE			

169RF

Knochenfräser, Stahl rostfrei
bone cutter, stainless steel

Handstück | Straight Handpiece

FIGURE	SHANK	REF	Ø
169RF	HPXL	330 106 412 297	023
		L mm	35,0
			2
APPLICATION & HYGIENE			

186RF

Vorkörner - Stahl rostfrei
initial bur - stainless steel

Handstück | Straight Handpiece

FIGURE	SHANK	REF	Ø
186RF	HP	330 104 684 377	018
Winkelstück Right Angle			
186RF	RA	330 204 684 377	018
		L mm	12,0
			5
APPLICATION & HYGIENE			

187RF

Vorkörner mit Stop - Stahl rostfrei
initial bur with Stop - stainless steel

Winkelstück | Right Angle

FIGURE	SHANK	REF	Ø
187RF	RA	330 204 685 377	018
		L mm	19,0
			5
APPLICATION & HYGIENE			

203RF

Spiralbohrer - Stahl rostfrei
twist drill - stainless steel

Winkelstück | Right Angle


FIGURE	SHANK	REF	Ø	
203RF	RAL	330 205 417 364	006	008
		L mm	7,0	7,0
			2	2
APPLICATION & HYGIENE				

203S

Spiralbohrer mit Stop zum Dekortikalisieren - Stahl rostfrei
twist drill with stop for decorticating - stainless steel

Winkelstück | Right Angle

FIGURE	SHANK	REF	Ø
203S	RA	330 204 449 336	012
		L mm	3,0
			2
APPLICATION & HYGIENE			

RCB00

Kieferkammkonturierer aus rostfreiem Stahl
ridge contouring bur of stainless steel

Winkelstück | Right Angle


FIGURE	SHANK	Ø
RCB00	RAL	075
		L mm
APPLICATION & HYGIENE		

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232RF

Scheibe, Umfang schneidend
disc, peripheral cutting

*in combination with mandrel 305.HP.050 Page 35

Unmontiert | Unmounted



FIGURE	SHANK	REF	Ø	
232RF	UM	330 900 320 101	190	220
		L mm	0,1	0,1
			10	10

PILOTBOHRER

PILOT DRILLS

1001

Winkelstück | Right Angle

6 mm
8 mm
10 mm
12 mm



FIGURE	SHANK	Ø					
A1001	RAXL	010					
B1001	RAXL		013				
C1001	RAXL			015			
D1001	RAXL				018		
E1001	RAXL					020	
F1001	RAXL						025
		L mm	15,0	15,0	15,0	15,0	15,0
			1	1	1	1	1
		APPLICATION & HYGIENE					
		ROTATION SPEED	200 - 600 rpm				

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
Winkelstück | Right Angle


FIGURE	SHANK	REF	Ø				
229XLRF	RAL	330 205 556 001	040	045	050	060	070
		L mm	18,0	18,0	18,0	18,0	18,0
			5,0	5,5	6,0	7,0	8,0
			4,0	4,5	5,0	6,0	7,0
			1	1	1	1	1
APPLICATION & HYGIENE							

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
Turbine | Friction Grip

FIGURE	SHANK	REF	Ø
GT48L	FG	700 314 287 484	016
		L mm	5,0
			1
APPLICATION & HYGIENE		   	

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Turbine | Friction Grip

FIGURE	SHANK	REF	Ø
GT135	FG	700 314 161 484	016
		L mm	8,0
			1
APPLICATION & HYGIENE		   	

HALS-, NASEN-, OHREN-CHIRURGIE INSTRUMENTE

EAR, NOSE, THROAT SURGERY INSTRUMENTS

4301

Bohrer rund, Stahl rostfrei
burs round, stainless steel

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BS48 FG-RA

24 x FG | 24 x RA

Dimensions (LxHxW)	81x51x31 mm
APPLICATION & HYGIENE	

BS4 RA

4 x RAL

autoklavierbar | autoclavable

Dimensions (LxHxW)	56.6x39x25.5 mm
APPLICATION & HYGIENE	

BS23 HP

23 x HP

autoklavierbar | autoclavable

Dimensions (LxHxW)	100x68x39,5 mm
APPLICATION & HYGIENE	

* Materialeignung der eingesetzten Mittel beachten (pH-Wert neutral, Eignung für Al)
 Observe suitability of the material for means used (pH value, suitability for Al)
 Tenir compte de l'adéquation des matériaux des produits utilisés (valeur pH neutre, adéquation pour Al)
 Es preciso tener en cuenta la idoneidad para el material del producto utilizado (valor pH neutro, adecuado para Al)

EDELSTAHL - BOHRERSTÄNDER, ROSTFREI
STAINLESS STEEL BUR BLOCKS, NON-CORROSIVE

BS12S FG-RA

12 x FG/RA

Dimensions (LxHxW)	60x40x45 mm
APPLICATION & HYGIENE	**

BS12S HP

12 x HP

Dimensions (LxHxW)	60x65x45 mm
APPLICATION & HYGIENE	**

BS24S HP

24 x HP

Dimensions (LxHxW)	120x65x45 mm
APPLICATION & HYGIENE	**

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INSTRUCTIONS FOR PROCESING & GENERAL APPLICATION AND SAFETY INSTRUCTIONS

Instructions for the processing (Cleaning, disinfection, and sterilization) of instruments from Jota AG Issued: January 2022

The medical devices produced and sold by Jota AG are re-usable unless their label contains explicit information to the contrary. However, as a rule, it is the sole responsibility of the doctor/expert using the devices to decide whether, depending on the respective case and the potential wear and tear of the products, he can re-use the products and how frequently he uses them. In case of doubt, it is always advisable to discard the products early and to replace them. The manufacturer Jota AG cannot guarantee the faultless function and performance of the products combined with a maximum degree of safety if the products are overused. These reprocessing instructions apply in principle to all medical devices making up the product range supplied by Jota AG. Any particular features and/or exclusions that may concern individual items or groups of items are referred to separately.

Fundamental points

All instruments are classified as semi-critical instruments and are to be cleaned, disinfected, and sterilized prior to each application; this is required as well for the first use after delivery of the unsterile instruments (cleaning and disinfection after removal of the protective packaging, sterilization after packaging). An effective cleaning and disinfection are indispensable requirements for an effective sterilization of the instruments.

You are responsible for the sterility of the instruments. Therefore, please ensure that only sufficiently device and product specifically validated procedures will be used for cleaning, disinfection, and sterilization, that the used devices (WD, sterilizer) will be maintained and checked regularly, as well as that the validated parameters will be applied for each cycle.

Please pay attention to avoid a higher contamination of the complete bur block during application; otherwise, it is necessary to clean and disinfect the bur block as well as all instruments inside (after removal).

Additionally, please pay attention to the legal provisions valid for your country as well as to the hygienic instructions of the doctor's practice or of the hospital. This applies particularly to the different guidelines regarding the inactivation of prions (not relevant for USA). It also refers to country specific specification of workplaces for cleaning and disinfection.

Some instruments require additional aspects. For this, pay attention-on to chapter "Specific aspects".

Cleaning and disinfecting

Basic:

If possible, an automated procedure (WD (Washer-Disinfector)) should be used for cleaning and disinfection of the instruments. A manual procedure – even in case of application of an ultrasonic bath – should only be used if an automated procedure is not available; in this case, the significantly lower efficiency and reproducibility of a manual procedure has to be considered. The pre-treatment step is to be performed in both cases.

Pre-treatment:

Please remove coarse impurities of the instruments directly after application (within a maximum of 1 h).

Procedure:

1. Rinse the instruments at least 1 min under running water (temperature < 35 °C/95 °F).
2. Soak the instruments at least for the given soaking time in the pre-cleaning solution (e.g. ID 212 forte, 2% for 5 Minutes; by the use of an ultrasonic bath) so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments. Assist cleaning by careful brushing with a soft brush (at least three times after beginning of soaking, aids see chapter "Specific aspects").
3. Then activate ultrasonic treatment for an additional soaking time (but not less than 5 min).
4. Then, remove the instruments of the pre-cleaning solution and post-rinse them at least three times intensively (at least 1 min) with cold tap water.
5. In case of still visible contamination repeat steps 2, 3, and 4, otherwise discard the instrument. This is especially relevant for diamond instruments.

Pay attention to following points during selection of cleaning detergent 1):

- > fundamental suitability for the cleaning of instruments made of metallic or plastic material
- > suitability of the cleaning detergent for ultrasonic cleaning (no foam development)
- > compatibility of the cleaning detergent with the instruments (see chapter „material resistance,,)

Pay attention to the instructions of the detergent manufacturer regarding concentration, temperature and soaking time as well as post-rinsing. Please use only freshly prepared solutions as well as only demineralized sterile or low contaminated water (max. 10 germs/ml) as well as low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water, and a soft, clean, and lint-free cloth and/or filtered air for drying, respectively.

- 1) In case of application of a cleaning and disinfection detergent for this (e.g. in consequence of personnel's safety) please consider, that this should be aldehyde-free (otherwise fixation of blood impurities), possess a fundamentally approved efficiency (for example VAH/DGHM or FDA/ EPA approval/clearance/ registration or CE marking), be suitable for the disinfection of instruments made of metallic or plastic material, and be compatible with the instruments (see chapter „material resistance,,). Please consider, that a disinfectant used in the pre-treatment step serves only the personnel's safety, but cannot replace the disinfection step later to be performed after cleaning.

Automated cleaning/disinfection (recommended; WD (Washer-Disinfector)):

Pay attention to following points during selection of the WD:

- > fundamentally approved efficiency of the WD (for example CE marking according to EN ISO 15883 or DGHM or FDA approval/clearance/ registration)
- > possibility for an approved program for thermal disinfection (AO value > 3000 or – in case of older devices - at least 5 min at 90°C/194 °F; in case of chemical disinfection danger of remnants of the disinfectant on the instruments)
- > fundamental suitability of the program for instruments as well as sufficient rinsing steps in the program
- > post-rinsing only with demineralized sterile or low contaminated water (max. 10 germs/ml, max. 0.25 endotoxin units/ml), for example purified/highly purified water

- > only use of filtered air (oil-free, low contamination with microorganisms and particles) for drying
- > regularly maintenance and check/calibration of the WD

Pay attention to following points during selection of the cleaning detergent:

- > fundamental suitability for the cleaning of instruments made of metallic or plastic material
- > additional application – in case of non-application of a thermal disinfection – of a suitable disinfectant with approved efficiency (for example VAH/DGHM or FDA/EPA approval clearance/registration or CE marking) compatible to the used cleaning detergent
- > compatibility of the used detergents with the instruments (see chapter „material resistance,,)

Pay attention to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing.

Procedure:

1. Transfer the pre-cleaned instruments in the WD by the use of a small piece's basket or tray.
2. Start the program:
3 min pre-rinsing with cold tap water (<30°C) drain
10 min cleaning at 50 - 60°C with 0.2-1% Neodisher® Mediclean Dental drain
1 min rinsing with demin. water (40-45°C) drain
1 min rinsing with cold demin. water (<30°C) drain
5 min thermal disinfection for 5 min at 90-92°C with demin. water drain
30 min automatic drying at 100°C (wd program)
3. Remove the instruments of the WD after end of the program.
4. Check and pack the instruments immediately after the removal (see chapters „check,, „maintenance,, and „packaging,, if necessary, after additional post-drying at a clean place).

The fundamental suitability of the instruments for an effective automated cleaning and disinfection was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of the WD G 7835 CD, Miele & Cie. GmbH & Co., Gütersloh, (thermal disinfection) and the pre-cleaning and cleaning detergent neodisher® Mediclean Dental (Dr. Weigert GmbH & Co. KG, Hamburg) considering to the specified procedure.

Manual cleaning and disinfection (not recommended):

Pay attention to following points during selection of the cleaning and disinfection detergents:

- > fundamental suitability for the cleaning and disinfection of instruments made of metallic or plastic material
- > in case of application of an ultrasonic bath: suitability of the cleaning detergent for ultrasonic cleaning (no foam development)
- > application of a disinfectant with approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) compatible with the used cleaning detergent
- > compatibility of the used detergents with the instruments (see chapter „material resistance,,)

Combined cleaning/disinfection detergents should not be used.

Only in case of extremely low contamination (no visible impurities) combined cleaning/disinfection could be used.

Pay attention to the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing. Please use only freshly prepared solutions as well as only demineralized sterile or low contaminated water (max. 10 germs/ml) as well as low endotoxin contaminated water (max. 0.25 endotoxin units/ml), for example purified/highly purified water, and a soft, clean, and lint-free cloth and/or filtered air for drying, respectively.

Manual Procedure after Pre-treatment:

Manual Cleaning

- Soak the instruments for the given soaking time in the cleaning solution (e.g. ID 212 forte, 2% for 5 Minutes; by the use of a ultrasonic bath) so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments. Assist cleaning by careful brushing with a soft brush (at least three times after at beginning of soaking, aids see chapter "Specific aspects").
- Then activate ultrasonic treatment for an additional soaking time (but not less than 15 min).
- Then, remove the instruments of the cleaning solution and post-rinse them at least three times intensively (at least 1 min) with water.
- Check the instruments (see chapters „check, and „maintenance“).

Manual Disinfection

- Soak the instruments for the given soaking time in the disinfectant solution (e.g. ID 212 forte, 2% for 1 Minutes; by the use of a ultrasonic bath) so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments.
- Then, remove the instruments of the disinfectant solution and post-rinse them at least five times intensively (at least 1 min) with water.
- Dry and pack the instruments immediately after the removal (see chapter „packaging,“ if necessary, after additional post-drying at a clean place).

The fundamental suitability of the instruments for an effective cleaning and disinfection was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of the pre- cleaning and cleaning detergent and the disinfectant ID 212 forte (Dürr Dental SE, Bietigheim-Bissingen) considering the specified procedure.

Check

Check all instruments after cleaning or cleaning/disinfection, respectively, on corrosion, damaged surfaces, and impurities. Do not further use damaged instruments (for limitation of the numbers of re-use cycles see chapter „reusability,“). Still dirty instruments or instruments with visible residuals (under regular light conditions) are to be cleaned and disinfected again.

Maintenance

Instrument oils or grease must not be use except for steel instruments. In that case use only instrument oils (white oil) admitted to steam sterilization considering the maximum possible sterilization temperature, with approved biocompatibility and without mono-, di, or triethanolamine as corrosion inhibitor.

Packaging

Please insert the cleaned and disinfected instruments in the corresponding bur blocks (if required) and pack them in single-use sterilization packaging (single or double packaging), which fulfill the following requirements (material/process):

- EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance)
- suitable for steam sterilization (temperature resistance up to at least 142 °C (288 °F), sufficient

steam permeability)

- sufficient protection of the instruments as well as of the sterilization packaging to mechanical damage

Sterilization

Please use for sterilization only the listed sterilization procedures; other sterilization procedures must not be applied.

Steam sterilization

- fractionated vacuum/dynamic air removal procedure^{2,3} (with sufficient product drying⁵)
- steam sterilizer according to EN 13060/EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- validated according to EN ISO 17665 (valid IQ/OQ (commissioning) and product specific performance qualification (PQ))
- maximum sterilization temperature 138 °C (280 °F; plus tolerance according to EN ISO 17665)
- sterilization time (exposure time at the sterilization temperature):

Area	fractionated vacuum/dynamic air removal	gravity displacement
USA	at least 4 min at 132 °C (270 °F), drying time at least 20 min ⁴	not recommended
other countries	at least 3 min ⁵ at 132 °C (270 °F) / 134 °C (273 °F), drying time at least 20 min ⁴	not recommended

² at least three vacuum steps

³ The less effective gravity displacement procedure must not be used in case of availability of the fractionated vacuum procedure, will require significantly longer sterilization times and is to be validated dependent on product, packaging, sterilizer, program, and parameters under sole responsibility of the user.

⁴ The effectively required drying time depends directly on parameters in sole responsibility of the user (load configuration and density, sterilizer conditions, ...) and by this is to be determined by the user. Nevertheless, drying times less than 20 min must not be applied.

⁵ respectively 18 min (inactivation of prions, not relevant for USA)

The fundamental suitability of the instruments for an effective steam sterilization was demonstrated by an independent, governmentally accredited and recognized (§ 15 (5) MPG) test laboratory by application of the steam sterilizer EHS 3870 (Tuttnauer, Breda Netherlands) and the fractionated vacuum/dynamic air removal procedure. For this, typical conditions in clinic and doctor's practice as well as the specified procedure were considered.

The flash/immediate use sterilization procedure must not be used. Do not use dry heat sterilization, radiation sterilization, formaldehyde, and ethylene oxide sterilization, as well as plasma sterilization.

Storage

Please store the instruments after sterilization in the sterilization packaging's at a dry and dust-free place.

Material resistance

Please take care that the listed substances are not ingredients of the cleaning or disinfection detergent:

- organic, mineral, and oxidizing acids (minimum admitted pH-value 5.5)
- strong lyes (maximum admitted pH-value 11, neutral/enzymatic or alkaline cleaner recommended)⁶
- organic solvents (for example: acetone, ether, alcohol, benzine)
- oxidizing agents (for example: hydrogen peroxide)
- halogens (chlorine, iodine, bromine)
- aromatic, halogenated hydrocarbons

⁶ For the bur blocks alkaline cleaners must not applied (maximum admitted pH-value 9).

Please do not clean any instruments and bur blocks by use of metal brushes or steel wool.

Please do not expose any instruments and bur blocks to temperatures higher than 142 °C (288 °F)!

Please do not apply acidic neutralizing agents or cleaning aids.

Reusability

The instruments can be reused – in case of adequate care and if they are undamaged and clean as indicated in chapter "Specific aspects". The user is responsible for each further use as well as for the use of damaged and dirty instruments (no liability in case of disregard).

Attachment A: Specific aspects

Diamond products and ceramic grinding tools:

- Use particular care during the cleaning of the grinding surfaces and ensure that all residues are removed. Repeat procedure, if visible residuals are given

Bur blocks/instrument trays:

- Cleaning and disinfecting only without products being loaded (products must not be cleaned and disinfected whilst they are in the bur block/instrument tray)

i

DE	Das Dokument in der jeweils aktuellsten und gültigen Version sowie in anderen Sprachen finden Sie auf unserer Webseite www.jota.ch .
EN	The document in other languages as well as the latest and applicable version you find on the website www.jota.ch .
FR	Vous trouverez les documents dans d'autres langues ainsi que la dernière version sur notre site www.jota.ch .
ES	Documentos en otras idiomas, así como también la versión más reciente encontrase en nuestra página web www.jota.ch .

Instrument group	brush	specific/additional procedure in case of						maximum admitted cycle number (confirmed by validation, but dependent on specific application)	recommended classification according to KRINKO/RKI/BfArM guidance (only German, with respect to intended use)
		pretreatment	manual cleaning/ disinfection	automated cleaning/ disinfection	maintenance	packing	sterilization		
stainless steel instruments	standard	standard	standard	standard	lubrication <u>not</u> admitted	standard	standard	10	critical B
regular steel instruments	standard	standard	standard	standard	lubrication <u>recommended</u>	standard	standard	10	critical B
silicone polisher	standard	standard	standard	standard	lubrication <u>not</u> admitted	standard	standard	5	critical B
endodontic instruments without stopper	endodontic brush	standard	standard	standard	lubrication <u>not</u> admitted	use of bur blocks not admitted	use of bur blocks not admitted	10	critical B
endodontic instruments with stopper	endodontic brush	mounted	mounted move the stopper at least three times during disinfection	mounted	lubrication <u>not</u> admitted	use of bur blocks not admitted	use of bur blocks not admitted	10	critical B
all other instruments	standard	standard	standard	standard	lubrication <u>not</u> admitted	standard	standard	10	critical B

General application and safety instructions for the medical device from Jota AG
Issued: January 2017

- Jota AG products (dental, maxillary surgery, general surgery,) must only be used by dentists, doctors and/or the respective experts who, due to their training and experience, are intensely familiar with the use of these products and who have the corresponding expertise in the respective specialist fields. The use of surgical products requires relevant expertise and experience in dental implantology, maxillary surgery and/or other surgical fields including diagnosis, preoperative planning and surgical techniques.
- It is the sole responsibility of the doctor in charge who, depending on the respective situation (indication), decides on the actual use of the Jota AG products for each individual case
- All JOTA AG products have been developed for specific applications. Therefore, inappropriate use can result in the premature wear and tear of the products and put patients and users at risk.

Application

- In order to avoid damaging the instruments, they must be removed from the blister pack by pulling off the back-sheet.
- It is essential to only use turbines as well as hand and angle pieces that are technically and hygienically faultless, maintained and cleaned.
- The instruments must be rotating when applied on material. They should not be placed on material and then brought to rotation.
- Rotating instruments need to be clamped as far down as possible with their speed set before applying them on the object. are used with the rotary instruments.
- Using the instruments for canting or leveraging should be avoided as it increases the risk of breakage.
- Depending on the application, it is recommended to use protective goggles while using the instruments. Users of diamond disks should use a disk protection device.
- Inappropriate use of the products leads to badly executed work and increased risk.
- When working with dry materials, it is recommended to use a suction cleaning device.
- In particular, users of hand tools should take care to use them gently and with consideration.
- The user must at all times avoid touching the instruments and parts without protection (protective gloves should be worn).
- Thermal bone damage caused by rotating and oscillating tools (e.g. pilot burr, conical burr, expansion burr) must at all times be avoided (user training, working at low speed and with sufficient cooling).
- During intraoral application attention has to be made to the fact that the products are protected against aspiration or falling on the floor.

Use of pressure

- Users of the instruments should at all times avoid applying excessive pressure. This can damage the working part of the instruments and cause the cutting edges to break off. At the same time, it generates excessive heat.
- The use of excessive pressure when using grinding tools can cause the abrasive particles to break off or the instrument to become clogged and lead to heat generation.
- During polishing, excess pressure can lead to heat generation.
- Due to overheating, excess pressure can damage the dental pulp or, due to broken off cutting edges, it can result in undesired rough surfaces. In such cases, even instrument breakage cannot be excluded.

Cooling

- In order to avoid excessive heat generation during preparation, a sterile water/sodium chloride solution supplied via a permanent external feeding device should be used to ensure sufficient cooling during use of the instruments.
- When using FG instruments that are more than 22 mm long or whose head diameter exceeds 2 mm, additional external cooling is required.
- Insufficient cooling will lead to irreversible damage to the bone and/or the adjacent tissue.

Storage, disinfection, cleaning and sterilization

- Unless there is explicit information to the contrary, all Jota AG products are supplied in non-sterile packaging and, depending on the application, they need to be sterilised prior to use. Prior to their first use on the patient and immediately after each use, all products need to be disinfected and sterilised. Inappropriate cleaning and sterilising of the instruments can result in the patient being infected with harmful bacteria.
- You will find detailed instructions for the disinfecting, cleaning and sterilising of products in the Instructions for the processing of instruments produced by Jota AG on the previous page. We would also be happy to provide you with these instructions at your request. They are also available on the internet at www.jota.ch.
- The products should be stored in appropriate, hygienically maintained containers. The same applies to sterilised instruments. The stored products must be protected from dust, humidity and recontamination. Instructions as to maximum storage duration must be adhered to.

Speed recommendations for rotary instruments

- Following the instrument-specific speed recommendations produces the best results.
- Exceeding the maximum admissible speed (rpm) when using long and pointed instruments tends to produce vibrations that can lead to the destruction of the instrument.
- When using working parts with diameters exceeding the thickness of the shaft, excessive speed can release great centrifugal forces that may cause the shaft to bend and/or the instrument to break. Therefore, the maximum admissible rpm must never be exceeded.

- Please consult the manufacturer's information (see catalogue or www.jota.ch) for the recommended and the maximum admissible speed ranges. Non-compliance with the maximum admissible speed puts safety at risk.
- Generally, the following rules apply:
 - The larger the working part of an instrument the lower the speed
 - Surgical instruments: suitable for geared down micro-motor hand and angle pieces 10:1 with stable ball bearings. Speed 600 to 800 rpm with physical and, possibly, sterile external cooling or internal cooling when using the respective hand piece.

Discarding worn instruments and parts

- Jota AG products can principally be reused several times – unless specifically indicated and labelled otherwise. Rotating instruments are subject to wear. The option of and accountability for multiple use of a product and the frequency of application is solely the decision and own responsibility of the treating clinician based on the application in each case and the possible wear of the products. If in doubt, the products should always be sorted out early and replaced.
- Broken off cutting edges of instruments cause vibrations and great forces of pressure, which, in turn, leads to broken preparation corners and rough surfaces.
- Bare patches on diamond instruments indicate a lack of abrasive particles and can be a sign of blunt instruments. This leads to excessive temperatures during instrument use.
- Instruments that are bent and/or do not run true should be discarded forthwith.
- With the reuse of disposable products the risk of infection cannot be excluded and a risk-free functional safety cannot be guaranteed.

Additional instructions regarding the use of trepans

- When using trepans, you have to proceed with particular care. For example, it is advisable not to exceed the recommended rpm speed ranges.
- In order to prepare for the actual use of a trepan, it should be set to produce counter-clockwise rotations creating a groove in the bone. Afterwards the trepan can be inserted into this groove and, using clockwise rotations, it can be moved further down.
- Carrying out a prior X-Ray is essential to establish the maximum possible drilling depth and to maintain the necessary distance, for example, to the mandibular nerve. As an additional safety measure to spare the nerve, the axial direction of the trepan countersink attachment, based on the sagittal level of the ascending branch, must be milled laterally at an angle of approx. 15-20°.

Further comments:

- Due to statutory regulations, returned goods can, on principle, only be accepted if the complete batch number is provided. This number can be found on the product packaging.

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OUR PERFORMANCE RECORDS | NOS RECORDS DE PERFORMANCE |
NUESTROS REGISTROS DE RENDIMIENTO

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Being one of the market leaders means constant further development in all corporate processes and compliance with all the regulatory requirements.

The laws and regulations relating to medical devices are among the most stringent worldwide, with the focus on patient protection.

High product quality is essential. Professionalism results from a high level of appreciation of the applications are from the qualification of the products.

- › Development and development verification
- › Clinical evaluations and biocompatibility evidence
- › Validation of processes
- › Batch management and proof obligations up to 30 years after the sale
- › Market surveillance

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All rights reserved. Reproduction, also by extract, are only permitted with written authorization of JOTA AG, Rüthi. All offers, orders and deliveries are subject to our „General Sales and Delivery Terms“. We reserve the right to modify our range of products and their design as well as to deviate from the illustrations and data shown.

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JOTA AG Rotary Instruments
Jota AG Rotary Instruments